

# UltraDose®

## GERMICIDAL Hospital Grade Ultrasonic Cleaner

Clean and disinfect in one easy step with dual action UltraDose® Hospital Grade Ultrasonic Cleaner Concentrate.

The unique blend of an ultrasonic cleaning solution and an E.P.A. registered hospital grade disinfectant come together in this effective product. UltraDose® Germicidal eliminates dangerous microorganisms while simultaneously ultrasonically cleaning and disinfecting instruments, apparatus and equipment. Help maintain the safety of your office and staff and keep your instruments clean and disinfected with UltraDose® Germicidal.

**Don't take infection control for granted!**

Available in 16 oz. containers.  
Each container makes  
16 gallons of solution.



## COVID-19 Approved\*

### \*Emerging Viral Pathogen Claim

UltraDose Germicidal Hospital Grade Ultrasonic Cleaner Concentrate has demonstrated effectiveness against viruses similar to COVID-19 Coronavirus on hard, non-porous surfaces. Therefore, UltraDose Germicidal Hospital Grade Ultrasonic Cleaner Concentrate can be used against COVID-19 Coronavirus when used in accordance with the directions for use against Human Coronavirus on hard, non-porous surfaces.

- E.P.A. registered Product
- For use in an ultrasonic cleaner
- Disinfects instruments and the ultrasonic bath
- Cleans non-porous instrument and appliance surfaces, apparatus and equipment
- Performs as an ultrasonic cleaning solution, disinfectant, deodorant, germicidal, fungicidal, virucidal and tuberculocidal all-in-one
- Concentrated; one ounce makes one gallon. Non-irritating and fragrance free



*because clean matters.*

# UltraDose®

## GERMICIDAL

### Hospital Grade Ultrasonic Cleaner



- ✓ Fragrance free
- ✓ Super concentrated
- ✓ Dual purpose proprietary formulation
- ✓ Promotes office staff safety

Use UltraDose® Germicidal to clean and disinfect instruments and appliances in an ultrasonic cleaning system and disinfect the ultrasonic bath at the same time.

UltraDose® Germicidal can also be used as a pre-soak to disinfect instruments prior to ultrasonic cleaning.

**Active Ingredient:** Phenylphenol


**Container Size:** 16 oz.

**Dilution:** 1 oz./gallon

**Gallons/container:** Makes 16 gallons of ready to use solution

**E.P.A. Registration #:** 34810-31-72869

**Product Code:** UD036



- ⊘ COVID-19\*
- ⊘ HIV-1
- ⊘ Pseudomonas Acrusinosa
- ⊘ Trichophyton
- ⊘ Salmonella
- ⊘ Staphylococcus
- ⊘ Seratia Marcescens
- ⊘ Herpes Simplex 2
- ⊘ Mycobacterium Bovis (TB)